

one step cleaner and neutral disinfectant

GS High Dilution Disinfectant 256 is a neutral disinfectant cleaner that offers economy and less frequent product container switch out. Just a ½ ounce of concentrate per gallon of water provides total cleaning, disinfection and deodorization in one labor saving, cost effective step.





Spartan Chemical Company, Inc. 1110 Spartan Drive Maumee, OH 43537 1-800-537-8990 www.spartanchemical.com

NON-ALKALINE - NO RINSING REQUIRED - LEAVES NO FILM

The neutral pH, quaternary disinfectant cleaner is mild, and formulated to kill a broad spectrum of microorganisms without rinsing or leaving a film residue.

MULTI-SURFACE USES

GS High Dilution Disinfectant 256 is versatile. It can clean and disinfect a variety of hard, nonporous, inanimate environmental surfaces including: floors, walls, metal surfaces, stainless steel surfaces, porcelain, glazed ceramic tile, plastic surfaces, bathrooms, shower stalls, bathtubs, and cabinets. GS High Dilution Disinfectant 256 is excellent for use in larger areas such as operating rooms and patient care facilities. This product effectively cleans floors without dulling the shine of floor finish.

GS High Dilution Disinfectant 256 can be used in a variety of businesses: healthcare, veterinary facilities, schools, child and adult day care, transportation terminals, office buildings, manufacturing facilities, lodging establishments, retail business, athletic/recreational facilities, sports stadiums, amphitheaters and convention centers.

KILLS GERMS IN MINUTES

Effective against the following Bacteria:

Pseudomonas aeruginosa PRD-10, Salmonella enterica subsp. enterica serotype choleraesuis, Staphylococcus aureus, Staphylococcus aureus (clinical isolate), Bordetella bronchiseptica, Corynebacterium ammoniagenes, Enterobacter aerogenes, Enterobacter cloacae, Enterobacter cloacae (clinical isolate), Enterococcus faecalis, Enterococcus faecalis, Enterococcus faecalis, Enterococcus faecalis (clinical isolate), Escherichia coli, Escherichia coli (clinical isolate), Fusobacterium necrophorum, Klebsiella pneumoniae subsp. pneumoniae, Lactobacillus

casei subsp. rhamnosus, Listeria monocytogenes, Pasteurella multocida, Proteus vulgaris, Proteus mirabilis ATCC 9921, Proteus mirabilis ATCC 25933, Salmonella enterica subsp. enterica serotype paratyphi B, Salmonella enterica subsp. enterica serotype typhi, Salmonella enterica subsp. enterica serotype typhimurium, Salmonella enterica subsp. enterica serotype pullorum, Serratia marcescens, Shigella sonnei, Shigella flexneri Type 2b, Shigella dysentariae, Staphylococcus aureus subsp. aureus, Staphylococcus epidermidis, Staphylococcus epidermidis (clinical isolate), Streptococcus pyogenes (Clinical - Flesh Eating Strain BIRD M3), Streptococcus pyogenes Group A, Xanthamonas maltophilia (clinical isolate), Vancomycin resistant Enterococcus faecalis (VRE), Methicillin resistant Staphylococcus aureus (MRSA), and Vancomycin intermediate resistant Staphylococcus aureus (VISA), Community Associated Methicillin Resistant Staphyloccus aureus (CA-MRSA) NRS 123 Genotype USA 400, Community Associated Methicillin Resistant Staphyloccus aureus (CA-MRSA) NRS 384 Genotype USA 300

Effective against the following Viruses:

HIV-1, HIV-2, Hepatitis B virus(HBV), Hepatitis C virus (HCV), Herpes Simplex Type 1 (causative agent of fever blisters), Herpes Simplex Type 2(genital), Influenza A2/Hong Kong, Vaccinia, Rotavirus, Human Coronavirus (ATCC VR-740, Strain 229E), SARS Associated Coronavirus, Rabies Virus, Bovine Viral Diarrhea Virus (BVDV), Pseudorabies, Bovine Rhinotracheitis, Feline Leukemia, Feline Picornavirus, Canine Distemper Virus, Avian Influenza A Virus, and Porcine Respiratory & Reproductive Syndrome Virus (PRRSV), Pandemic 2009 H1N1 Influenza A virus

Fungicidal claims:

Candida albicans, Trichophyton mentagrophytes

DIRECTIONS FOR USE

Use GS High Dilution Disinfectant 256 for these application and more:

GENERAL CLEANING: Apply product to soiled area with a mop, cloth, sponge, or hand pump trigger sprayer. Product can be diluted at ½ ounce per gallon or applied at full strength. Thoroughly wet soiled surface and scrub as necessary.

DISINFECTION: To disinfect inanimate, hard, non-porous surfaces add ½ ounce of this product per gallon of water. Apply solution with a mop, cloth, sponge, or hand-pump trigger sprayer so as to wet all surfaces thoroughly. Allow to remain wet for 10 minutes, then remove excess liquid.

For heavily soiled areas, a pre-cleaning step is required. Prepare a fresh solution for each use.

See product label for complete directions for use and more applications.

Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting www.spartanchemical.com.

SPECIFICATIONS

pH (Concentrate): 6.5-7.5

Color: Clear

Scent: Fragrance-Free

Stability

• Shelf @ 24°C/75°F: One year

• Freeze/Thaw: 3 cycles

Contains biodegradable surfactants. Phosphate Free. EPA Reg. No. 1839-167-5741





Label copy is provided in English and Spanish. Secondary labels are also available.

GUARANTEE: Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture. Use product as directed and read all precautionary statements.