



Plexiglass Care

EQUIPMENT



Microfiber Cloth

KEY

- 1 Glass Cleaner
- 2 Disinfectant/Sanitizer
- 3 All Purpose Cleaner

RECOMMENDED PRODUCTS

Spartan recommends these daily cleaners and disinfectants.



1 BioRenewables® Glass Cleaner RTU (323903)



2 Hepacide® Quat II (325103)



2 Sani-Tyze® (319503)



3 Sign Off® (329003)



3 Xcelente® Multi-Purpose Cleaner (603100)



3 Xcelente® Multi-Purpose Cleaner RTU Handi Spray® (318603)

CLEANING AND DISINFECTING INFECTION BARRIERS (SNEEZE GUARD)

- 1 **CLEAN:** Using a microfiber cloth, wipe down plexiglass with the cleaner of your choice as needed. Most glass and multi-surface cleaners and wipes are suitable for plexiglass with the exception or high pH or ammoniated products.
- 3
- 2 **DISINFECT:** When selecting a disinfectant or disinfecting wipe, look for a neutral pH (6.5 - 7.5). Noting that these products may leave a residue behind and require an additional wipe with a clean, dry microfiber cloth after disinfectant dwell time is achieved.

Always follow guidance from your plexiglass manufacturer on cleaning recommendations.

